


# CARDIOPLEGIA SYSTEM



**Indication:** Indicated for preconditioning and application of cardioplegic solutions (blood or crystalloid) to the heart, for a safe and protected cardiac stoppage during an extracorporeal procedure.

### TECHNICAL SPECIFICATIONS - Heat Exchanger

Raw Material	Polycarbonate
Capillary Material	Polyester
Max Flow	500ml/min
Effective Heat Exchange Area	0.14m <sup>2</sup>
Prime Volume	44ml
Max Water Temperature	42°C
Max Water Pressure	40psi
Solution Inflow/Outflow Connector	3/16"
Water Inflow/Outflow Connector	Hansen
Air Extraction Connector	Luer Lock

MODEL	RESERVOIR	CODE
Crystalloid	RC 2000	209649
Incor	RC 500	220288
4:1	n/a	220172
Adult Cardiomicroplegia	n/a	220154

4:1



Crystalloid



Incor



Cardiomicroplegia

